< relaxmyway

QUATTROMED III

Limited Edition

< relaxmyway

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< relaxmyway

Since the company was founded in the year 2000, the name CASADA has stood for high grade products in the area of wellness and fitness.

Today, CASADA – products are available worldwide. What distinguishes CASADA products is their continuous development, unmistakable design and their absolute focus on improving quality of life as well as their reliable and excellent quality.

Congratulations!

With the purchase of this massage device you have shown awareness of your health.

In order for you to enjoy the advantages of your massage device in the long term, we ask you to read and adhere to the instruction manual carefully.

We wish you a lot of fun with your personal Quattromed III Limited Edition.

Table of Contents

12	Safety instructions
12	Start up
13	Equipment / Functions
13	Troubleshooting
14	Specials
16	Remote Control
16	Technical Data
17	Storage / Care
17	Conformity declaration
17	Warranty conditions

Safety Instructions

Please carefully read the following instructions before using this massage device to ensure its proper function and optimal effect. Please retain these instructions for use!

- This massage device complies with the recognised principles of technology and current safety regulations.
- This equipment is maintenance free. Any necessary repairs must be performed by an authorised professional. Improper use and unauthorised repairs are prohibited for safety reasons and will void the warranty.
- Never touch the power plug with wet hands.
- Please avoid the equipment coming into contact with water, high temperatures and direct sunlight.
- Be sure to unplug the equipment from the power point during extended periods of non-use to eliminate the risk of a short circuit.
- Do not use damaged plugs, cords or loose sockets. Damaged plug or cords must be replaced by the manufacturer, service representative or qualified personnel. In the event of malfunction unplug the equipment immediately. We assume no liability for damages resulting from wrongful or improper use.
- In order to prevent excessive stress on muscles and nerves, the recommended massage time of 30 minutes should not be exceeded.
- Never use pointy or sharp objects.
- Any massage – also massage by hand – has to be avoided during pregnancy and/or in case of one or more complaints in the area of the massage: recent injuries, thrombotic illnesses, inflammations and swelling of any kind as well as cancer. Consultation of your physician is recommended prior to any treatment of afflictions and ailments.
- Do not use the device while driving.
- If you depend on electrical aids such as a pacemaker, please consult your physician prior to usage in any case.

Start Up

Carefully loosen the screw on the back of the massage device. Remove the cover and connect the plug to the adapter. Attention! The screw is not to be used again. Just remove it.

Please connect the adapter with the device and subsequently plug the power plug into the power socket. Now turn on the device by the remote control.

Equipment / Functions

In the new Quattromed III Limited Edition the massage heads normally made from plastic are substituted with heads of jade. The natural powers of jade are increased by an integrated carbon heating system. The massage will be intensified by sending healing warmth into the depths of the muscles.

The heat is permeating through the stone slowly and is only noticeable after a few minutes. During this heating phase the tissue is already being relaxed by the rotating movement of the massage heads and surface hardening and tensions are being released. The optional vibration also allows for the stimulation of the circulation and lymphatic system. Blood circulation is effectively improved and the purification facilitated.

Please position the Quattromed III Limited Edition on a stable chair. The burden of your bodyweight intensifies the massage effect. Please ensure that your whole body is positioned centrally so that the massage heads run parallel to your spine. Using the rolling function you can vary the width between the massage heads, pressing the according button.

The massage function will switch off automatically after 15 minutes. In order to increase the lifespan of the motors, please switch on the massage heads again only after a minimum of 30 minutes have passed. This gives the motors time to cool down again. You should also allow a period of rest for your body. We recommend conducting the massages only two, maximum three times a day in order to not overstress your muscles.

For your own safety we equip the machines with a protection against overheating. Heat is created in the motors when placed under great stress and excessive pressure. Before this could create any risk, the machines will switch off automatically and will only allow switching on again after an appropriate downtime between 30 and 40 minutes.

Troubleshooting

Error

The massage device cannot operate even though it has been connected to the power supply.

Adjustment

- Check the connection of the adapter.
- The device could be overheated after prolonged use. It will automatically cool down to resume its normal condition. Please turn off the power and allow the device to cool off for 30 to 50 minutes.

Specials

Specials

KNEADING

Skin and muscles are either clasped and kneaded between the thumb and the index finger or with both hands. This massage technique is used in particular to ease tension.

ROLLING

Pleasant alternation between relaxation and releasing the muscles – also perfectly suited for a final massage, for relaxation and recovery of the muscles.

CARBON HEATING

Deeply penetrating carbon infrared heat, blood vessels expand, the blood flow increases and the muscle tension is reduced. In addition, it carries out a highly beneficial effect on nerve cells, which are responsible for reporting pain to the brain.

VIBRATION

The vibration setting moves parts of the body with equal vibrations, which penetrate deeply into the muscles. This stimulates the circulation in the muscles.

REVERSE

The reverse function changes the direction of the massage heads.

JADE HEATING

Jade is a natural store and conductor of heat. These high-quality jade stones, purposefully integrated in the massage heads dispense deep warmth and a relaxing effect.

AUTOMATIC

Specifically developed automatic programmes for the whole body, which offer a big range of different massage techniques. Additionally you can also turn on the heating function, if you'd like.

CAR ADAPTER

The handy car adapter makes your massage device portable.

SPOT

The massage heads can massage a particular spot.

WIDTH


The distance between the both massage heads can be set (narrow, medium and wide).

Remote Control

- [1] On/Off
Switch the device on/off here. When switched off, the massage stops and the mechanism returns automatically to the lower rest position.
- [2] Carbon and Jade Heating
On/off additional carbon and jade heating.
- [3] Time (5 – 15 minutes)
The preset massage time is 15 minutes.
This button sets the desired massage duration.
(5 minute intervals / 5 – 10 – 15).

SELECT MESSAGE AREA

- [4] Switch the massage on/off over the entire back area.
- [5] Switch the massage on/off in the shoulder and neck area.
- [6] Switch the massage on/off in the back and lumbar area.
- [7] Vibration (in seat area)
Start/stop the vibration massage.
- [8] Kneading
This button switches the kneading function on/off.
- [9] Reverse (only possible in kneading mode)
This button determines the rotation direction of the massage heads.
- [10] Rolling
This button switches the rolling function on/off.
- [11] Width (only possible in rolling mode)
This button changes the spacing of the massage heads during a rolling massage.
- [12] Spot (only possible in kneading mode)
When held down, these buttons move the position of the massage heads upwards and downwards.
When the respective button is released, the mechanism stops at the desired point.

Technical Data Quattromed III LE

Measurements:	72 x 20 x 45 cm
Weight:	5,6 kg
Power source:	Input: AC 100 -240 V ~ 50/60 Hz, 1600 mA Output: DC 12 V, 5000 mA max. 60 Watt
Rated power:	max. 60 Watt
Runtime:	15 minutes
Certificates:	

Storage / Care

Regular care will extend the life of your Quattromed III Limited Edition whilst also ensuring safe operation at all times.

1. Recommended storage: dry, free of dust.
2. Keep liquids and open fire away from the Quattromed III Limited Edition
3. Clean the equipment using a dry, soft cloth.
4. Always unplug the Quattromed III Limited Edition when not in use.

Conformity Declaration

Corresponds to the European guidelines EMC 2014/30/EU and LVD 2014/35/EU.

Warranty Conditions

The dealer grants the statutory warranty period on the product described overleaf. The warranty applies from the date of purchase. The date of purchase is demonstrated by the sales receipt.

The dealer is to repair or replace products registered within Germany which are found to be defective free of charge. Parts subject to wear and tear such as cloth covers are not covered. To assert your claim, simply return the defective product with the sales receipt before the end of the warranty period.

The warranty claim does not apply if it is established that a defect has occurred e.g. through external influence or as a result of repair or alteration not carried out by the manufacturer or an authorised dealer.

The warranty granted by the seller is limited to repair or replacement of the product. Under this warranty, the manufacturer or seller has no further liability and is not responsible for damage caused by disregard of the operating instructions and/or improper use of the product.

If inspection of the product by the seller reveals that the warranty claim relates to a defect not covered by the warranty or that the warranty period has expired, the cost of the inspection and repair are to be borne by the customer.

Legal note: All contents do not contain healing statements. The diagnosis and treatment of illnesses and other physical disorders require a treatment by a physician, non-medical practitioner or therapist. The statements are strictly informative and may not be used as a substitute for a professional medical treatment. Every user is beholden to ascertain by accurate examination and, if required, after consultation of a specialist, if the application is the right one for their specific cases. Each dose or application takes place at the risk of the user. We dissociate ourselves from any healing statements or promises.

Casada International GmbH
Obermeiers Feld 3
33104 Paderborn, Germany
info@casada.com
www.casada.com

Copyright Images and texts are subject to the copyright of Casada International GmbH and may not be used for other purposes without express confirmation.
Copyright © 2017 Casada International GmbH. All rights reserved.

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Copyright © 2017 CASADA International GmbH.
Les images et textes publiés sont soumis à la protection des droits d'auteur, ils ne peuvent être extraits sans autorisation formelle. CASADA International GmbH est le titulaire des droits d'auteur. Copyright © 2017 CASADA International GmbH.